

BEAUTY

CATALOGUE



PRODUCING QUALITY GOODS FOR 15 YEARS

We are Producers of quality surgical, dental, beauty and general stainless steel products in singleuse and multiuse type. We guarantee on our products and facilitate our customers with easy return policies.



For more information visit www.mjwintl.net



Manufacturers and Exporters of Surgical, Dental and Beauty Instruments

Head Office

BI24S 88/8 Bara Pathar
Christian Town, Sialkot
Pakistan (51310)
0092 524 264829

Branch Office

Flat 1, Kinsale Grange
57 Langley Park Road,
Sutton-Surrey SM2-5GZ UK.
Mob: 00447701084895

Contact Details

Email: info@mjwintl.net
Website: www.mjwintl.net
Mobile: 0092 321 712 1510
Mobile: 0092 333 371 3443



Nail
nippers

Instruments for easy shortening and forming strong finger & toe nails.

The different sizes and cutting head forms demonstrate the broad range of manicure or pedicure applications.

There are two different joints available. Lap joint nippers have the nipper shanks lying one on top of the other at the joint. Box joint nippers have one shank put through the other nippers shank.

While one can cut and form the finger nails according to desire and fashion, it is recommended to cut and form toe nails relatively short and straight. Toe nails should not reach over the end of the toes.



MJW-301-110
Nail nipper
lap joint
wire spring
chequered handle
Size 10cm



MJW-302-110
Nail nipper
lap joint
wire spring
plain handle
Size 10cm



MJW-303-110
Nail nipper
lap joint
wire spring
fluted handle
Size 10cm



MJW-310-120
Nail nipper
lap joint
wire spring
Chequered handle
Size 12cm



MJW-311-120
Nail nipper
lap joint
wire spring
plain handle
Size 12cm



MJW-320-120
 Corner Cutter
 box joint
 fine concave jaws
 single spring with lock
 Size 12cm



MJW-325-130
 Corner Cutter
 box joint
 normal jaws
 single spring
 Size 12cm



MJW-328-130
 Corner Cutter
 box joint
 fine jaws
 single spring
 Size 13cm



MJW-330-130
 Corner Cutter
 box joint
 fine/concave jaws
 single spring
 Size 13cm



MJW-335-130
 Corner Cutter
 box joint
 single spring
 extra fine point
 Size 13cm



MJW-336-130
 Corner Cutter
 box joint
 single spring
 extra fine point
 Size 13cm



MJW-340-135
 Corner Cutter
 box joint
 single spring with lock
 Size 13.5cm



MJW-340-135
 Corner Cutter
 box joint
 double spring
 delicate jaws
 with one cutter blade
 Size 13.5cm