

SURGICAL

CATALOGUE



PRODUCING QUALITY GOODS FOR 15 YEARS

We are Producers of quality surgical, dental, beauty and general stainless steel products in singleuse and multiuse type. We guarantee on our products and facilitate our customers with easy return policies.



For more information visit www.mjwintl.net

Manufacturers and Exporters of Surgical, Dental & Beauty Instruments

Head Office

BI24S 88/8 Bara Pathar
Christian Town, Sialkot
Pakistan (51310)
0092 524 264829

Branch Office

Flat 1, Kinsale Grange
57 Langley Park Road,
Sutton-Surrey SM2-5GZ UK.
Mob: 00447701084895

Contact Details

Email: info@mjwintl.net
mjwinternational2005@gmail.com
Mobile: 0092 321 712 1510
Mobile: 0092 333 371 3443



FACTORY OVERVIEW



INTRODUCTION

for fifteen years MJW Surgical has been supplying the Healthcare Industry with finest quality surgical instruments & innovative products for surgery. With our extensive knowledge and experience in the medical field, our products are specifically designed to offer both exceptional quality and value for money to Healthcare Professionals and Organisations. During the last few years there have been numerous developments within our company, with the Introduction of many new innovative products and services to our portfolio, which has allowed MJW Surgical to make a significant impact on the healthcare market. Our head office is based in Sialkot, Pakistan housing all Customer Service and Administrative functions. Additionally, we have an office in the UK covering Europe and the USA. Our range of specialised surgical instrumentation is uniquely manufactured in Sialkot to the highest standards of craftsmanship. Our goal is to continue to develop new, innovative medical products that will enhance the quality of surgical procedures.

OUR SERVICE COMMITMENT

QUALITY of All surgical instruments and products manufactured and distributed by MJW Surgical conform to the highest European and International Standards with all products bearing the relevant CE mark for their classification.

GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any product proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts.

QUOTATION SERVICE

We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.

PRODUCT INFORMATION

Sectional catalogues are available for most surgical specialities.

INNOVATION

Our objective is to work in close conjunction with all medical and nursing personnel to develop new and innovative products and services.

29-112
For aspirating fluid from the excretory ducts of large glands and for keeping the operating area clear 170 mm

29-113
Trocar
For puncturing
150 mm

29-102

	A	B	C
29-102	5 mm	1 mm	4.0 mm
29-103	6 mm	1 mm	5.0 mm
29-103A	7 mm	1 mm	5.5 mm
29-103B	8 mm	1 mm	6.0 mm
29-103C	9 mm	1 mm	6.5 mm
29-103D	10 mm	2 mm	7.5 mm
29-103E	12 mm	2 mm	10 mm
29-103F	14 mm	2 mm	11 mm

Uterine suction tube

29-104
29 cm
Specimen & aspirating curette

29-105
29 cm
4.5 mm
Specimen & aspirating curette

Aspirating currettes
29 cm

29-106
8 mm

29-107
8 mm

29-108
10 mm

29-109
12 mm

29-110
14 mm

29-111
12 cm
Coupling piece

29-115 150mm
Suction handle, with vacuum regulator and finger cut-off, brass, mirror-finish chromium-plated

29-116 170mm
Suction handle, with hand regulation valve and quick-shutter, brass, mirror-finish chromium-plated

olive
Ø 9 mm

olive
Ø 9 mm

Both suction handles have bayonet fitting and can be used for the suction cannulae.

29-114
For clearing large quantities of fluid
190 mm

29-117
Trocar
For puncturing
170 mm

29-118
Hollow needle,
sharp
150 mm

29-119
double curved
For general surgery
230 mm


29-120
For tonsillectomy,
oral and laryngeal
operations
190 mm

29-121
For tonsillectomy,
oral and laryngeal
operations
220 mm

Ø 11 mm

Ø 7.5 mm

With bayonet fitting to be used with suction handles (29-115 - 29-116)




29-122
For tonsillectomy,
oral and laryngeal
operations
180 mm

29-123
Suction cannula
for tonsillectomy,
oral and laryngeal
operations
220 mm

29-124
Suction cannula
for tonsillectomy,
oral and laryngeal
operations
180 mm

29-125
S-shaped blunt cannula
for aspirating fluid from
the excretory ducts of
large glands and for
keeping the operating
area clear
170 mm



29-126
POOLE
Suction cannula,
with outer tube
to be screwed off
225 mm

29-127
POOLE
Suction cannula,
with outer tube
to be screwed off
225 mm



29-128
YANKAUER
with especially large
central bore, 4.8 mm
diam.
285 mm

29-129
YANKAUER
with central bore,
2.0 mm diam.
285 mm

29-130
DE-BAKEY
with central bore,
3.0 mm diam. with
lateral slots
275 mm

29-131
DE-BAKEY
with finger cut-off
and stylet
270 mm



29-132
DE-BAKEY
with finger cut-off
and stylet
160 mm

29-133
with finger
cut-off
190 mm



29-134 - 137
PLESTER
with finger cut-off and stylet, thumb piece parallel to curve, with LUER hub
195 mm

29-138 - 142
FERGUSSON
with finger cut-off and stylet, thumb piece 90° to curve, with LUER hub
180 mm

29-143 4,0 mm
FERGUSSON
insulated, with finger cut-off and stylet, thumb piece 90° to curve, with LUER hub
180 mm

29-144 3,0 mm
29-145 5,0 mm
FRAZIER
with finger cut-off and stylet, thumb piece 90° to curve, with LUER hub
180 mm

29-146 5,0 mm diam. Tip 6,5 mm wide, Nasal suction tube, with LUER hub
180 mm

29-147 4,5 mm
29-148 5,5 mm
Nasal suction tubes
175 mm

29-149
GORNEY
Suction dissector
180 mm

29-150 4 mm wide
29-151 5 mm wide
Suction raspatories, sharp, bayonet-shaped, with stylet
200 mm

29-152
LOPEZ-REINKE
Tonsil dissector with suction device
235 mm

29-153
WULLSTEIN
Suction tube, conical, silver, olive with LUER hub
160 mm

29-154 1,5 mm-29-155 2,0 mm
WULLSTEIN
Suction tube, cylindrical, silver, with stylet olive with LUER hub
160 mm

YASARGIL usable length diam.		YASARGIL usable length diam.	
29-156 130 mm	1,5 mm	29-164 130 mm	2,5 mm
29-157 150 mm	1,5 mm	29-165 150 mm	2,5 mm
29-158 180 mm	1,5 mm	29-166 180 mm	2,5 mm
29-159 220 mm	1,5 mm	29-167 220 mm	2,5 mm
29-160 130 mm	2,0 mm	29-168 130 mm	3,5 mm
29-161 150 mm	2,0 mm	29-169 150 mm	3,5 mm
29-162 180 mm	2,0 mm	29-170 180 mm	3,5 mm
29-163 220 mm	2,0 mm	29-171 220 mm	3,5 mm
		29-172 180 mm	4,5 mm



29-173 - 175
LANDOLT
150 MM
diam 3 mm

29-173 8 mm / curved
29-174 11mm / curved
29-175 15mm / curved

29-176 500 mm-29-177 350 mm
with rubber tip, with LUER hub

29-178 2,0 mm diam
with finger cut-off, thumb piece 90° to curve, with LUER hub
250 mm

29-179 2,5 mm diam
29-180 3,0 mm diam



TUNGESTEN CARBIDE INSTRUMENTS
 FOR THOSE.....
 WHO DEMAND SECURITY & ACCURACY.



29-210
 GENERAL BIOLOGY DISSECTING SET
 SINGLE LEATHER CASE.



29-211
 BASIC DISSECTING SETS
 VINYLE POUCH.



29-212
 ANATOMY DISSECTING SET
 DOUBLE FOLDING LEATHER CASE.



29-213
 LARGE ANIMAL SET
 DOUBLE FOLDING
 LEATHER CASE.



29-214
 MICRO DISSECTING SET
 DOUBLE FOLDING
 LEATHER CASE.



29-215
 NEW INSTRUCTORS DISSECTING SET
 (Set Of 12)



29-216
 BASIC DISSECTING SETS
 In Attractive Zipper Case
 (set of 9)



29-217
 BASIC DISSECTING SETS (Set of 13)
 (In Attractive Zipper Case)



Tungsten Carbide is One Of The Hardest Alloys. We Bond These Inserts Into High Quality German Stianless SteelVery Special Reasonns Of The End Users. We Offer A Wide Range Of Instruments Which Includes Scissors, Forceps And Needle Holders Of All kinds Made Of German Stainless Steel. Briefly, We Strive To End Up With Quality, Precision And Conformity To Our Buyers/specification And Standard By All Means. The Process Indeed Involves Complex Machining Operations, Skill And Experience. Let Us Please Secure Your Atteion To Serve You Better. **MANAGEMENT**