

SURGICAL

CATALOGUE



PRODUCING QUALITY GOODS FOR 15 YEARS

We are Producers of quality surgical, dental, beauty and general stainless steel products in singleuse and multiuse type. We guarantee on our products and facilitate our customers with easy return policies.



For more information visit www.mjwintl.net



Manufacturers and Exporters of Surgical, Dental & Beauty Instruments

Head Office

BI24S 88/8 Bara Pathar
Christian Town, Sialkot
Pakistan (51310)
0092 524 264829

Branch Office


Flat 1, Kinsale Grange
57 Langley Park Road,
Sutton-Surrey SM2-5GZ UK.
Mob: 00447701084895

Contact Details

Email: info@mjwintl.net
mjwinternational2005@gmail.com
Mobile: 0092 321 712 1510
Mobile: 0092 333 371 3443



FACTORY OVERVIEW



INTRODUCTION

for fifteen years MJW Surgical has been supplying the Healthcare Industry with finest quality surgical instruments & innovative products for surgery. With our extensive knowledge and experience in the medical field, our products are specifically designed to offer both exceptional quality and value for money to Healthcare Professionals and Organisations. During the last few years there have been numerous developments within our company, with the Introduction of many new innovative products and services to our portfolio, which has allowed MJW Surgical to make a significant impact on the healthcare market. Our head office is based in Sialkot, Pakistan housing all Customer Service and Administrative functions. Additionally, we have an office in the UK covering Europe and the USA. Our range of specialised surgical instrumentation is uniquely manufactured in Sialkot to the highest standards of craftsmanship. Our goal is to continue to develop new, innovative medical products that will enhance the quality of surgical procedures.

OUR SERVICE COMMITMENT

QUALITY of All surgical instruments and products manufactured and distributed by MJW Surgical conform to the highest European and International Standards with all products bearing the relevant CE mark for their classification.

GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any product proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts.

QUOTATION SERVICE

We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.

PRODUCT INFORMATION

Sectional catalogues are available for most surgical specialities.

INNOVATION

Our objective is to work in close conjunction with all medical and nursing personnel to develop new and innovative products and services.

TUNGESTEN CARBIDE INSTRUMENTS
FOR THOSE.....
WHO DEMAND SECURITY & ACCURACY.



29-210
GENERAL BIOLOGY DISSECTING SET
SINGLE LEATHER CASE.



29-211
BASIC DISSECTING SETS
VINYLE POUCH.



29-212
ANATOMY DISSECTING SET
DOUBLE FOLDING LEATHER CASE.



29-213
LARGE ANIMAL SET
DOUBLE FOLDING
LEATHER CASE.



29-214
MICRO DISSECTING SET
DOUBLE FOLDING
LEATHER CASE.



29-215
NEW INSTRUCTORS DISSECTING SET
(Set Of 12)



29-216
BASIC DISSECTING SETS
In Attractive Zipper Case
(set of 9)



29-217
BASIC DISSECTING SETS (Set of 13)
(In Attractive Zipper Case)




Tungsten Carbide is One Of The Hardest Alloys. We Bond These Inserts Into High Quality German Stianless SteelVery Special Reasonns Of The End Users. We Offer A Wide Range Of Instruments Which Includes Scissors, Forceps And Needle Holders Of All kinds Made Of German Stainless Steel. Briefly, We Strive To End Up With Quality, Precision And Conformity To Our Buyers/specification And Standard By All Means. The Process Indeed Involves Complex Machining Operations, Skill And Experience. Let Us Please Secure Your Atteion To Serve You Better. **MANAGEMENT**

Instruments with Tungsten Carbide Inserts

TC GOLD



STANDARD
30-101 to 30-106




30-101 14.5 cm
30-102 16 cm
30-103 18 cm
Straight

30-104
30-105
30-106
Curved



STANDARD
30-107 to 30-112




30-107 14.5 cm
30-108 16 cm
30-109 18 cm
Straight

30-110
30-111
30-112
Curved



STANDARD
30-113 to 30-118




30-113 14.5 cm
30-114 16 cm
30-115 18 cm
Straight

30-116
30-117
30-118
Curved



MAYO
30-119 to 30-124



30-119 14.5 cm
30-120 17 cm
30-121 23 cm
Straight

30-122
30-123
30-124
Curved



MAYO-LEXER
30-125 to 30-126
16 cm



30-125
30-126
Straight
Curved



METZENBAUM
30-127 to 30-133



30-127 14.5 cm
30-128 18 cm
30-129 20 cm
30-130 23 cm
30-131 25 cm
30-132 28 cm
30-133 30 cm
Curved



METZENBAUM
30-134 to 30-137



30-134 10 cm
30-135 12 cm
Straight

30-136 10 cm
30-137 12 cm
Curved



IRIS
30-138 to 30-139




30-138 Straight
30-139 Curved
11.5 cm

Instruments with Tungsten Carbide Inserts

TC GOLD



MAYO-STILLE
30-150 to 30-153




30-150 15 cm
30-151 17 cm
Straight

30-152
30-153
Curved



METZENBAUM-FINO
30-140 to 30-149



30-140 14 cm
30-141 16 cm
30-142 18 cm
30-143 20 cm
30-144 23 cm
Straight

30-145
30-146
30-147
30-148
30-149
Curved



LISTER
30-154 to 30-159

30-154 9 cm
30-155 11 cm
30-156 14 cm
30-157 16 cm
30-158 18 cm
30-159 20 cm



BRUNS
30-160
24 cm Serrated


Instruments with Tungsten Carbide Inserts
TC GOLD




Instruments with Tungsten Carbide Inserts
TC GOLD



Instruments with Tungsten Carbide Inserts
TC GOLD



Instruments with Tungsten Carbide Inserts
TC GOLD




Instruments with Tungsten Carbide Inserts
TC GOLD



90-456
18 cm / 7"
For soft wires up to 2 mmØ
For hard wires up to 1.5 mmØ



90-457
22 cm / 8 1/2"
For soft wires up to 3.5 mm
For hard wires up to 2.5 mm



TC wire cutter For wire
90-463 38.5 cm / 15" 5 mm
90-464 45 cm / 17 1/2" 6 mm
90-465 56 cm / 22" 7.5 mm



90-458
22 cm / 8 1/2"
For soft wires up to 2.5 mmØ
For hard wires up to 2.2 mmØ



90-459 14 cm / 5 1/2"
90-460 16 cm / 6 1/2"
90-461 18 cm / 7"
90-462 22 cm / 8 1/2"
For soft wires up to 3 mmØ
For hard wires up to 2.2 mmØ



SECTION PAGE	
Surgical Holloware	31
	414 to 424