

DENTAL

CATALOGUE



PRODUCING QUALITY GOODS FOR 15 YEARS

We are Producers of quality surgical, dental, beauty and general stainless steel products in singleuse and multiuse types. We guarantee on our products and facilitate our customers with easy return policies.



For more information visit www.mjwintl.net



Manufacturer and Exporter of Surgical, Dental and Beauty Instruments

Head Office

BI24S 88/8 Bara Pathar
Christian Town, Sialkot
Pakistan (51310)
0092 524 264829

Branch Office


Flat 1, Kinsale Grange
57 Langley Park Road,
Sutton-Surrey SM2-5GZ UK.
Mob: 00447701084895

Contact Details

Email: info@mjwintl.net
mjwinternational2005@gmail.com
Mobile: 0092 321 712 1510
Mobile: 0092 333 371 3443



FACTORY OVERVIEW



INTRODUCTION

for fifteen years MJW Surgical has been supplying the Healthcare Industry with finest quality surgical instruments & innovative products for surgery. With our extensive knowledge and experience in the medical field, our products are specifically designed to offer both exceptional quality and value for money to Healthcare Professionals and Organisations. During the last few years there have been numerous developments within our company, with the introduction of many new innovative products and services to our portfolio, which has allowed MJW Surgical to make a significant impact on the healthcare market. Our head office is based in Sialkot, Pakistan housing all Customer Service and Administrative functions. Additionally, we have an office in the UK covering Europe and the USA. Our range of specialised surgical instrumentation is uniquely manufactured in Sialkot to the highest standards of craftsmanship. Our goal is to continue to develop new, innovative medical products that will enhance the quality of surgical procedures.

OUR SERVICE COMMITMENT

QUALITY of All surgical instruments and products manufactured and distributed by MJW Surgical conform to the highest European and International Standards with all products bearing the relevant CE mark for their classification.

GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any product proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts.

QUOTATION SERVICE


We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.

PRODUCT INFORMATION

Sectional catalogues are available for most surgical specialities.

INNOVATION


Our objective is to work in close conjunction with all medical and nursing personnel to develop new and innovative products and services.



Stainless Steel Instrument Trays
Made of 18/8 Stainless Steel



ALUMINIUM INSTRUMENT TRAYS



Dental Instrument Tray
786-1855
150x100x10 mm
185x105x10 mm
200x100x10 mm
200x150x15 mm
280x190x20 mm
290x190x10 mm

Scalar Tray
786-1856
190x100x15 mm
Plain Base
190x100x15 mm
Grooved Base

Instrument Tray
786-1857
200x100x50 mm
200x150x50 mm
210x115x50 mm
220x150x50 mm
240x180x50 mm
250x200x50 mm
300x150x50 mm
300x200x50 mm
300x250x50 mm
please contact for more sizes

Instrument Tray with Lid
786-1858
200x100x50 mm
200x150x50 mm
210x115x50 mm
220x150x50 mm
240x180x50 mm
250x200x50 mm
300x150x50 mm
300x200x50 mm
300x250x50 mm
please contact for more sizes

Instrument Box with Lid
786-1859
170x70x30 mm
175x75x40 mm
180x60x30 mm
180x80x40 mm
180x90x40 mm
200x400x40 mm
200x100x50 mm
220x120x50 mm
230x150x50 mm
250x125x50 mm
260x150x50 mm
300x150x50 mm
Available with & without locks
please contact for more sizes

Tray With Lid
786-1860
200x75x35 mm
210x115x45 mm

Kidney Tray 786-1861 150 mm 200 mm 250 mm 300 mm	Sterilizing Drum 786-1863 100x100 mm 120x120 mm 150x150 mm 160x120 mm 190x140 mm 200x200 mm 240x120 mm 240x160 mm 240x240 mm	Gallipot 786-1864 2 OZ 4 OZ 6 OZ 8 OZ
Lotion Bowl 786-1862 80x40 mm 100x50 mm 120x55 mm 150x55 mm 180x60 mm		Forceps Jar 786-1870 50x100 mm 50x170 mm



All the above impression trays are available in these 5 sizes

UPPER	U1	U2	U3	U4	U5
A	62	60	53	54	52
B	83	75	73	69	68

All Dimension are in mm

LOWER	L1	L2	L3	L4	L5
A	63	61	57	60	53
B	82	81	78	74	71

All Dimension are in mm



786-1881/10
Perforated Partial Set of 10

786-1891/8
Solid Full Denture Set of 8



All the above impression trays are available in these 4 sizes

UPPER		XL	L	M	S
	A	70	67	64	60
UPPER		XL	L	M	S
	B	79	76	70	73

All Dimension are in mm

LOWER		XL	L	M	S
	A	64	64	60	57
LOWER		XL	L	M	S
	B	86	70	73	70

All Dimension are in mm



786-1899/8
Perforated Full Denture Set of 8

786-1907/6
Solid Set of 6



All the above impression trays are available in these 3 sizes

UPPER		U1	U2	U3
	A	67	64	57
UPPER		U1	U2	U3
	B	81	71	70

All Dimension are in mm

LOWER		L1	L2	L3
	A	60	57	54
LOWER		L1	L2	L3
	B	81	79	76

All Dimension are in mm




786-1913/6
Perforated Set of 6

786-1919/8
Solid Full Denture Set of 8



All the above impression trays are available in these 4 sizes



UPPER		XL	L	M	S
	A	70	67	64	60
	B	79	76	70	73

All Dimension are in mm

LOWER		XL	L	M	S
	A	64	64	60	57
	B	86	70	73	70

All Dimension are in mm



786-1927/8
Perforated Full Denture Set of 8

English Pattern

786-1935/8
Set of 8

Mesh type




786-1936/8
Set of 8
Punched type



786-1937/8
Set of 8
Solid type



All the above impression trays are available in these 4 sizes




UPPER	Size	A	B	C
	0	67	50	12
	1	72	51	13
	2	78	55	15
	3	80	58	15

All Dimension are in mm

LOWER	Size	A	B	C
	0	70	50	14
	1	77	51	15
	2	79	55	15
	3	80	60	17

All Dimension are in mm



Impression trays for crown & bridge work



786-1938
U1P

Solid

786-1939
L1P



786-1940
U1P

Perforated

786-1941
L1P

English Pattern
786-1942/8
Solid Full Denture
Set of 8



786-1942
U0

786-1943
U1

786-1944
U2

786-1945
U3



786-1946
L0

786-1947
L1

786-1948
L2

786-1949
L3

Upper Sizes 0,1,2,3,
Depth of the tray 15,9 mm
Height of palate 12,5 mm

Lower Sizes 0,1,2,3,
Depth of the tray 12,9 mm



786-1950
U0

786-1951
U1

786-1952
U2

786-1953
U3



786-1954
L0

786-1955
L1

786-1956
L2

786-1957
L3

786-1950/8
Perforated Full Denture Set of 8

Stainless Steel Impression Trays
Abdruckloeffel aus edelstahl
Cubetas para impresiones de acero inoxidable



786-1958
Adjustable
Perforated



786-1959
Adjustable
Solid



786-1960
For upper right jaw
For lower left jaw

Perforated



786-1961
For upper left jaw
For lower right jaw

Available in 4 Sizes
1.....44 mm
2.....55 mm
3.....60 mm
4.....65 mm

Also available in Solid pattern

LAB INSTRUMENTS

