

DENTAL

CATALOGUE



PRODUCING QUALITY GOODS FOR 15 YEARS

We are Producers of quality surgical, dental, beauty and general stainless steel products in singleuse and multiuse types. We guarantee on our products and facilitate our customers with easy return policies.



For more information visit www.mjwintl.net



Manufacturer and Exporter of Surgical, Dental and Beauty Instruments

Head Office

BI24S 88/8 Bara Pathar
Christian Town, Sialkot
Pakistan (51310)
0092 524 264829

Branch Office


Flat 1, Kinsale Grange
57 Langley Park Road,
Sutton-Surrey SM2-5GZ UK.
Mob: 00447701084895

Contact Details

Email: info@mjwintl.net
mjwinternational2005@gmail.com
Mobile: 0092 321 712 1510
Mobile: 0092 333 371 3443



FACTORY OVERVIEW



INTRODUCTION

for fifteen years MJW Surgical has been supplying the Healthcare Industry with finest quality surgical instruments & innovative products for surgery. With our extensive knowledge and experience in the medical field, our products are specifically designed to offer both exceptional quality and value for money to Healthcare Professionals and Organisations. During the last few years there have been numerous developments within our company, with the introduction of many new innovative products and services to our portfolio, which has allowed MJW Surgical to make a significant impact on the healthcare market. Our head office is based in Sialkot, Pakistan housing all Customer Service and Administrative functions. Additionally, we have an office in the UK covering Europe and the USA. Our range of specialised surgical instrumentation is uniquely manufactured in Sialkot to the highest standards of craftsmanship. Our goal is to continue to develop new, innovative medical products that will enhance the quality of surgical procedures.

OUR SERVICE COMMITMENT

QUALITY of All surgical instruments and products manufactured and distributed by MJW Surgical conform to the highest European and International Standards with all products bearing the relevant CE mark for their classification.

GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any product proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts.

QUOTATION SERVICE

We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.

PRODUCT INFORMATION

Sectional catalogues are available for most surgical specialities.

INNOVATION

Our objective is to work in close conjunction with all medical and nursing personnel to develop new and innovative products and services.



786-2033
Waldsachs
15 cm



786-2034
Universal
16 cm



786-2035
Nance
13.5 cm



786-2036
Aderer
14.5 cm



786-2037
Modell marburg
Serrated Jaws
14 cm



786-2038
Modell marburg
Smooth jaws
14 cm



786-2039
Modell marburg
Carbid Inserts jaws
14 cm



786-2040
Goslee
One Serrated Jaws
14 cm



786-2041
Aderer
12 cm



786-2042
Aderer- Heavy
12 cm



786-2043
Tweed
14.5 cm



786-2044
O'Brien
13.5 cm



786-2045
Wire bending pliers
12.5 cm



786-2048
Young
14 cm





786-2046
Adams
12.5 cm



786-2047
Frevert
12.5 cm



786-2049
Clasp forming
14 cm



**786-2075
How Plier**

Firmly hold the bracket tie wing for seating of the band. Carbide serrated tips for non slip grip. 49.50




**786-2076
De Bonding Plier**

Bracket / Adhesive remover with replaceable tool steel tip. Wedge between both edges of the base of the bracket and tooth surface. Easily lifting off the bracket or Bonding material. 49.50



**786-2077
Jara Back Plier**

A Plier for precise bending of wires. Precision grooves ensure accurate bending and closing of loops. Serrated tip for up to .020" wires. 49.50



**786-2078
Posterior Band Remover Long**

The padded tip rest on the Occlusal of the molar and the carbide tipped side removes the bracket or Bonding material. 49.50



**786-2079
Hard Wire Cutter**

Cuts hard wire up to .020" carbide insert for longer life. 49.50



**786-2080
Posterior Band Remover Small**

49.50



**786-2081
Pin & Ligature Cutter**

Standard size lips designed to cut soft wire up to .051". 49.50



**786-2082
Tweed Loop Forming Plier**

3 Size carbide tip allows to make 3 size loops. Can hold wire firmly upto .022" x .028". The other side serrated tip gives none slip grip. 49.50



All Pliers Also Available In Box Joint

**786-2083
Bird Beak Plier**

1/2" Long tip carbide inserted round tip and flat ground tip gives non slip grip for wire up to .030". 49.50



**786-2084
Angled Direct Bond Bracket Remover**

60° Angled tip allows for easy access to difficult anterior and posterior areas. 49.50



**786-2085
Safety hold, small Distal end Cutter**

Narrow and light. 49.50



**786-2086
Weingart Utility Plier**

Matching deep serrated pads holds wire firmly at a convex working angle. Tapered narrow tip fits easily between brackets. 49.50



**786-2087
Safety hold, Distal end Cutter**

Carbide tips for longer life and sharpness cuts and hold wire up to .022" x .028". 49.50




**786-2088
Three Prong Clasps Adjusting Plier**

Rounded jaws close and adjust clasps without niking wire upto .030". 49.50



**786-2089
Mini Pin & Ligature Cutter**

Narrow tips cuts pins, Ligatures and elastic smoothly and soft wire up to .015". 49.50



**786-2090
Light Wire Plier**

Wire forming Plier upto .061" the longer tapered beaks make it easier to bend small diameter loops. 3Grooves on square tip. 49.50



All Pliers Also Available In Box Joint



Rectangular Arch plier
13 cm / 5"
786-2091 Box Joint
786-2092 Screw Joint
786-2093 L-Key Joint
Short jaws, Tweed style, for use with square or rectangular wires. Blades parallel at 0.02", 0.51mm opening.



Weingart style plier
13 cm / 5"
786-2094 Box Joint
786-2095 Screw Joint
786-2096 L-Key Joint
Classical design, rounded tips. Useful for removal of ceramics.



Arch contouring plier
13 cm / 5"
786-2097 Box Joint
786-2098 Screw Joint
786-2099 L-Key Joint
for easily forming arches.



Arch forming & contouring plier
13 cm / 5"
786-2100 Box Joint
786-2101 Screw Joint
786-2102 L-Key Joint
With 3 grooves. For wire 0.016", 0.018" & 0.022".



Arch forming plier
13 cm / 5"
786-2103 Box Joint
786-2104 Screw Joint
786-2105 L-Key Joint
Without grooves



Micro Distal end cutter
13 cm / 5"
786-2106 Box Joint
786-2107 Screw Joint
786-2108 L-Key Joint
With safety hold, cuts arch wires. Cuts rectangular wire up to 0.56x0.71mm, round wire up to 0.51mm and twisted wire up to 0.41mm



Hard wire cutter
13 cm / 5"
786-2109 Box Joint
786-2110 Screw Joint
786-2111 L-Key Joint
For cutting orthodontic wires for fixed techniques any size. Cuts rectangular wire up to 0.56x0.71mm and round wire up to 0.51mm



Lindblom
13 cm / 5"
786-2121 Box Joint
786-2122 Screw Joint
786-2123 L-Key Joint



45° Ligature cutter
13 cm / 5"
786-2112 Box Joint
786-2113 Screw Joint
786-2114 L-Key Joint
45° angled for easy access to posterior areas. Cuts ligature wire up to 0.36mm



Miniature ligature cutter
13 cm / 5"
786-2115 Box Joint
786-2116 Screw Joint
786-2117 L-Key Joint
For soft wires. ligature wire up to 0.014" (0.36mm)

Distal wire cutter
15 cm / 5"
786-2118 Box Joint
786-2119 Screw Joint
786-2120 L-Key Joint



786-2124
Band seating & removing instruments

786-2125

786-2126

786-2127
Elastic Instruments

786-2128
Mershon Band Pusher
14 cm

786-2129
Band Setter with lin insert
15 cm

786-2130
13.5 cm

786-2131
14 cm

Bracket holding tweezers



786-2132
Separating pliers for separators
15.5 cm

ORTHODONTIC INSTRUMENTS



TUNGSTEN
CARBIDE INSTRUMENTS

