

SURGICAL

CATALOGUE



PRODUCING QUALITY GOODS FOR 15 YEARS

We are Producers of quality surgical, dental, beauty and general stainless steel products in single use and multi use type. We guarantee on our products and facilitate our customers with easy return policies.



For more information visit www.mjwintl.net



Manufacturers and Exporters of Surgical, Dental & Beauty Instruments

Head Office

BI24S 88/8 Bara Pathar
Christian Town, Sialkot
Pakistan (51310)
0092 524 264829

Branch Office

Flat 1, Kinsale Grange
57 Langley Park Road,
Sutton-Surrey SM2-5GZ UK.
Mob: 00447701084895

Contact Details

Email: info@mjwintl.net
mjwinternational2005@gmail.com
Mobile: 0092 321 712 1510
Mobile: 0092 333 371 3443



FACTORY OVERVIEW



INTRODUCTION

for fifteen years MJW Surgical has been supplying the Healthcare Industry with finest quality surgical instruments & innovative products for surgery. With our extensive knowledge and experience in the medical field, our products are specifically designed to offer both exceptional quality and value for money to Healthcare Professionals and Organisations. During the last few years there have been numerous developments within our company, with the introduction of many new innovative products and services to our portfolio, which has allowed MJW Surgical to make a significant impact on the healthcare market. Our head office is based in Sialkot, Pakistan housing all Customer Service and Administrative functions. Additionally, we have an office in the UK covering Europe and the USA. Our range of specialised surgical instrumentation is uniquely manufactured in Sialkot to the highest standards of craftsmanship. Our goal is to continue to develop new, innovative medical products that will enhance the quality of surgical procedures.

OUR SERVICE COMMITMENT

QUALITY of All surgical instruments and products manufactured and distributed by MJW Surgical conform to the highest European and International Standards with all products bearing the relevant CE mark for their classification.

GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any product proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts.

QUOTATION SERVICE

We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.

PRODUCT INFORMATION

Sectional catalogues are available for most surgical specialities.

INNOVATION

Our objective is to work in close conjunction with all medical and nursing personnel to develop new and innovative products and services.





14-110
IVES-FANSLER
Anuscope
145 mm
24mm Ø



14-111
KELLY
Sphincteroscope
160 mm
21mm Ø



14-112
KELLY
Sphincteroscope
160 mm
27mm Ø



14-113
KELLY
Proctoscope
140 mm
22mm Ø



14-114



14-115
YEOMAN




14-116
HERGET



14-117
Turnable
240 mm


14-118
Turnable
400 mm




Handle for (14-117, 14-118)






With Drawing Cut,
Attachments only


















15-110
MATHIEU
Foreign body forceps
Length of shaft
190 mm



15-111
MILLIN
Lobe grasping
forceps
235 mm



15-112
Ureter forceps
240 mm



15-113
MILLIN
Lobe grasping
forceps
245 mm



15-114
MILLIN
Bladder neck spreader
285 mm



15-115
GUYON
Catheter guide
(introducer)
365 mm



VEN BUREN
15-115-8 to 15-115-36
Set of 29

Charr. (French)
15-116 6
15-117 7
15-118 8
15-119 9
15-120 10
15-121 11
15-122 12
15-123 13
15-124 14
15-125 15
15-126 16
15-127 18
15-128 20

Charr. (French)
15-129 8
15-130 9
15-131 10
15-132 11
15-133 12
15-134 13
15-135 14
15-136 15
15-137 16
15-138 17
15-139 18
15-140 19
15-141 20
15-142 21
15-143 22
15-144 23
15-145 24
15-146 25
15-147 26
15-148 27
15-149 28
15-150 29
15-151 30

15-116 - 128
Female catheters,
metal
160 mm



15-129 - 151
DITTEL
Urethral bougies
and dilating sounds
345 mm



15-152 - 188
GUYON
Urethral bougies
and dilating sounds,
without filiform guides,
with thread in tip
255 mm



Fr.	Fr.	Fr.
15-189-08	15-189-19	15-189-30
15-189-09	15-189-20	15-189-31
15-189-10	15-189-21	15-189-32
15-189-11	15-189-22	15-189-33
15-189-12	15-189-23	15-189-34
15-189-13	15-189-24	15-189-35
15-189-14	15-189-25	15-189-36
15-189-15	15-189-26	
15-189-16	15-189-27	
15-189-17	15-189-28	
15-189-18	15-189-29	

Assisted by two haemostats or ligatures the foreskin is pushed over the bell-shaped protection - rod and the glass is placed into it.
By now, the wing - nut is tightened as close as possible, so the foreskin is in anular compressed. The foreskin is cut with a scalpel, right next to the compression ring. The circumcision clamp remains a few minutes to avoid bleeding in case that any large vein, specially the frenular artery, is cut through, a ligature must be set.
The advantage of the method described is the cosmetically favourable scar which is left after rejection of the small necrotized portion.



15-190 - 192	WINKELMANN
15-190	12mm
15-191	16mm
15-192	20mm



MOGEN
15-192A
7 cm



15-190 A	8 mm
15-190 B	11 mm
15-190 C	13 mm
15-190 D	16 mm
15-190 E	21 mm
15-190 F	26 mm
15-190 G	32 mm

Luer Lock



15-193
KLIMA ROSEGGER
35 X 1,8 mm
Sternal puncture needle



15-194
Universal trocar
with 4 tips

15-195 - 199 NELSON	Inside diam	Outside diam
15-195	2,2 mm	3,0 mm
15-196	3,0 mm	4,0 mm
15-197	4,0 mm	5,0 mm
15-198	5,0 mm	6,0 mm
15-199	5,5 mm	6,5 mm

15-200 - 203	Inside diam	Outside diam
15-200	2,2 mm	3,0 mm
15-201	3,0 mm	4,0 mm
15-202	4,0 mm	5,0 mm
15-203	5,0 mm	6,0 mm

for abdominal puncture,
with stop-cock and lateral
drain



15-204 - 205
BUELAU
Inside diam Outside diam
15-204 6,5 mm 8,0 mm
15-205 9,0 mm 10,0 mm
Empyema trocar



15-206 - 209
OCHSNER
Inside diam Outside diam
15-206 3,0 mm 4,0 mm
15-207 3,5 mm 4,5 mm
15-208 4,3 mm 5,3 mm
15-209 5,0 mm 6,0 mm
Gall bladder trocars

Record Cone 1

15-210
PIERCE
Antrum trocar
ø 4,0 mm



Luer Lock Cone

15-211
COAKLEY
Antrum trocar
ø 3,0 mm

Circumcision Clamps