

SURGICAL

CATALOGUE



PRODUCING QUALITY GOODS FOR 15 YEARS

We are Producers of quality surgical, dental, beauty and general stainless steel products in single use and multi use type. We guarantee on our products and facilitate our customers with easy return policies.



For more information visit www.mjwintl.net



Manufacturers and Exporters of Surgical, Dental & Beauty Instruments

Head Office

BI24S 88/8 Bara Pathar
Christian Town, Sialkot
Pakistan (51310)
0092 524 264829

Branch Office

Flat 1, Kinsale Grange
57 Langley Park Road,
Sutton-Surrey SM2-5GZ UK.
Mob: 00447701084895

Contact Details

Email: info@mjwintl.net
mjwinternational2005@gmail.com
Mobile: 0092 321 712 1510
Mobile: 0092 333 371 3443



FACTORY OVERVIEW



INTRODUCTION

for fifteen years MJW Surgical has been supplying the Healthcare Industry with finest quality surgical instruments & innovative products for surgery. With our extensive knowledge and experience in the medical field, our products are specifically designed to offer both exceptional quality and value for money to Healthcare Professionals and Organisations. During the last few years there have been numerous developments within our company, with the introduction of many new innovative products and services to our portfolio, which has allowed MJW Surgical to make a significant impact on the healthcare market. Our head office is based in Sialkot, Pakistan housing all Customer Service and Administrative functions. Additionally, we have an office in the UK covering Europe and the USA. Our range of specialised surgical instrumentation is uniquely manufactured in Sialkot to the highest standards of craftsmanship. Our goal is to continue to develop new, innovative medical products that will enhance the quality of surgical procedures.

OUR SERVICE COMMITMENT

QUALITY of All surgical instruments and products manufactured and distributed by MJW Surgical conform to the highest European and International Standards with all products bearing the relevant CE mark for their classification.

GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any product proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts.

QUOTATION SERVICE

We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.


PRODUCT INFORMATION

Sectional catalogues are available for most surgical specialities.

INNOVATION

Our objective is to work in close conjunction with all medical and nursing personnel to develop new and innovative products and services.

Orthopaedic Instruments 310 To 392		
		
SECTION	PAGE	
Bone Punches And Rongeurs	22	310-318
Amputation, Elevators, Raspatories, Bone Levers Bone Files And Curettes	23-24	319-347
Mallets, Tampers, Chisels, Gauges, Bone Holding, Rongeurs & Cutting Forceps, Meniscus Instruments	25-26	348-379
Wire Extension, Finger Nail instruments	27	380-392



90°

upwards cutting

85 mm

22-104
4 mm22-105
5 mm

1/4

22-104 - 105
KERRISON

90°

upwards cutting

140 mm

22-106
3,3 mm22-107
4,0 mm

1/4

22-106 - 107
HAJEK-KOFLER

130°

upwards cutting

KERRISON-JACOBY

22-108
1 mm

FERRIS-SMITH-KERRISON

22-109
2 mm22-110
3 mm22-111
4 mm22-112
5 mm

180 mm

1/4

1/4

22-108 - 112

90°

upwards cutting

KERRISON-JACOBY

22-113
1 mm22-114
2 mm

FERRIS-SMITH-KERRISON

22-115
3 mm22-116
4 mm22-117
5 mm

180 mm

1/4

1/4

22-113 - 117

Angled positions for rongeurs



22-118
3,0 mm
22-119
3,6 mm
22-120
4,2 mm
22-121
4,8 mm
22-122
5,8 mm



22-118 - 122
WEIL-BLAKESLEY



22-123
TAKAHASHI
3 mm

130°



22-124
TAKAHASHI
3 mm

LOVE-GRUENWALD

22-125
2 X 10 mm

22-126
3 X 10 mm

SPURLING

22-127
4 X 10 mm



22-125 - 127

22-128
2 mm

22-129
3 mm

22-130
4 mm

22-131
5 mm



22-128 - 131
CASPAR

150°

22-132
2 mm

22-133
3 mm

22-134
4 mm



22-132 - 134
CASPAR

150°
Lower part of jaws movable

22-135
2 mm

22-136
3 mm

22-137
4 mm



22-135 - 137
CASPAR

Length of shaft 140 mm	Length of shaft 160 mm	Length of shaft 185 mm	size of jaws
22-138	22-143	22-148	2 X 12 mm
22-139	22-144	22-149	3 X 12 mm
22-140	22-145	22-150	4 X 14 mm
22-141	22-146	22-151	5 X 14 mm
22-142	22-147	22-152	6 X 16 mm

150°

Length of shaft 140 mm	Length of shaft 160 mm	Length of shaft 185 mm	size of jaws
22-153	22-156	22-159	3 X 12 mm
22-154	22-157	22-160	4 X 14 mm
22-155	22-158	22-161	5 X 14 mm



22-153 - 161
CASPAR

150°
Lower part of jaw movable

Length of shaft 140 mm	Length of shaft 160 mm
22-162	22-164
22-163	22-165



22-162 - 167
CASPAR

22-168
2 X 10 mm
CUSHING

22-169
3 X 10 mm
LOVE-GRUENWALD

22-170
4 X 10 mm
SPURLING



150°

22-171
3 X 10 mm
LOVE-GRUENWALD

22-172
4 X 10 mm
SPURLING



22-171 - 172





ADSON
23-101
15 cm
Drill guide



23-102
14 cm



23-103
14 cm



23-104
30 cm
Straight ball spike



PERTHES
23-105
18 cm



PERTHES
23-105A
21 cm



BUTCHER
23-105B
Blade 45 cm
Post mortem saw



23-105C
Blade 40 cm
Post mortem saw



23-105D
Blade 20 cm
Post mortem saw



23-105A For children
23-105B For adults



23-106
PERCY
removable handles



23-107
PERCY
reversible handles



23-108
CHARRIERE
Length of blade
215 mm
complete with one saw
blade each(23-109 and 23-110)
350 mm

Spare saw blades, carbon steel
23-109 12 mm
23-110 20 mm



23-111
RUST
Length of blade
235 mm
complete with one saw
blade each(23-113, 23-114
and 2 saw blades 23-115)
390 mm
23-112 Bow only

Spare saw blades, carbon steel
23-113 6 mm
23-114 8 mm
23-115 12 mm
23-116 20 mm



23-117
BIER
Length of blade
270 mm
complete with one saw
blade each(23-118, 23-119 and 23-120)
420 mm

Spare saw blades, carbon steel
23-118 4 mm
23-119 6 mm
23-120 8 mm



23-121
CHARRIERE
Length of blade
180 mm
Spare saw blade,
carbon steel (23-122)
270 mm

23-123
CHARRIERE
Length of blade
240 mm
Spare saw blade,
carbon steel(23-124)
340 mm





23-125
SATTERLEE
Length of blade
200 mm
Spare saw blade,
carbon steel(23-126)
290 mm








23-127
LANGENBECK
Length of blade
115 mm
Metacarpal saw
235 mm






































can be sterilized in
a steam sterilizer upto
3 bar over pressure
(143 C).














26-101
Reposition forceps,
for small fragments
135 mm



26-102
HAASE
Reposition forceps,
for mandibular
fractures
165 mm



26-103
REILL
Reposition forceps,
170 mm



26-104
Reposition forceps
205 mm



26-105
MEYER
Reposition forceps
for fingers, with drill
guide for KIRSCHNER
wires up to diam. 1.0 mm
170 mm



26-106
Patella forceps
185 mm



26-107
Forceps with
ratchet fixation
140 mm



26-108
Forceps with
thread fixation
170 mm



26-109
Forceps with
thread fixation
235 mm



26-110
DINGMANN modif
185 mm



26-111
Sequestrum
forceps
200 mm

26-112
VAN BUREN
Sequestrum
forceps
225 mm


26-113
LANGENBECK
205 mm



26-114
FERGUSSON
205 mm



26-115
SEMB
190 mm



FARABEUF
26-116 19 cm
26-116A 23 cm
26-116B 26 cm
Adjustable

26-117
FARABEUF
255 mm



FARABEUF-LAMBOTTE
26-118 19 cm
26-118A 23 cm
26-118B 26 cm
Adjustable



26-119
KERN-LANE
240 mm




26-120
LANE
330 mm




26-121
KERN-LANE
240 mm



LANE
26-122
30 cm
26-122A
32.5 cm
Bone holding forceps.
Standard model, without ratchet.
26-122B
32.5 cm
Standard model, with ratchet.
26-122C
42.5 cm
With ratchet



HEY GROVES
26-123
42.5 cm
With ratchet
21 cm
26-123A
25 cm
26-123B
30.5 cm





26-138 - 139
Synovectomy rongeur
width of jaws 1.2 mm
130 mm



26-140 - 141
Synovectomy rongeur
width of jaws 2.0 mm
130 mm



26-142
STELLBRINK
Synovectomy rongeur
175 mm



26-146
LUER
150 mm



26-147
LUER
155 mm



26-148
BLUMENTHAL
155 mm



26-143
STELLBRINK
Synovectomy rongeur
170 mm



26-144
LUER-FRIEDMANN
145 mm



26-145
140 mm



26-149
BLUMENTHAL
150 mm



26-150 - 151
JANSEN
175 mm



26-152
BANE-HARTMANN
175 mm




26-153
ALEXANDER
140 mm



26-154
JANSEN
190 mm



26-155 - 156
LEMPERT
200 mm



26-160 - 161
LUER
175 mm



26-162
LUER
180 mm



26-163
LUER
180 mm



26-157
BEYER
180 mm



26-158
MEAD
165 mm



26-159
MEAD
165 mm



26-164
OLIVECRONA
205 mm



26-165 - 166
BOHLER
155 mm



26-167
MARQUARDT
200 mm



26-168
RUSKIN
190 mm



26-169
ROETTGEN-RUSKIN
240 mm



26-170 - 171
ECHLIN
230 mm





26-172
STILLE
235 mm



26-173
LEKSELL-STILLE
240mm





26-174
STILLE-RUSKIN
240mm







McIndoe
26-191A
17.5 cm



26-191
LISTON
curved on flat
240 mm



26-192
LISTON
straight
280 mm



26-193
LISTON
curved on flat
280 mm



SEMB
26-194A
24 cm



26-194
LISTON-STILLE
265 mm



26-195
LISTON-KEY-HORSLEY
(TUDOR EDWARDS)
255 mm





27-101 - 107
KIRSCHNER
Wire traction bows,
to be used with
KIRSCHNER-wires
up to 2.5 mm diam.,
With 3 hooks,
stainless steel



27-108
Wrench, double - ended,
for wire traction bows 27-101 - 27-115



27-109
Wire tightener, with
clamping jaws



27-110
Traction hook only,
large size



27-111
Traction hook only,
small size




Kuntscher
26 cm
27-111A 8 mm
27-111B 10 mm
27-111C 12.5 mm



Kuntscher
27-111D
35 cm




27-112
Fixation clamp, to secure clamping
jaws of 27-109 stainless steel



27-113
Wire traction bow, adjustable,
with 3 hooks, small size

27-114
Wire traction bow, adjustable,
with 3 hooks, large size



27-115
Wire traction bow, adjustable,
with 3 hooks, small size for
elbow and calcaneus, for
adolescents



27-116 - 118
KIRSCHNER
Wire traction bows, modified,
adjustable in width, holds wire
securely, with 1 hook



27-119
Fixation disc, slotted



27-120
idem, aluminum, transparent
to X-rays



27-121
Fixation disc, slotted



27-122
Screw driver, for above
fixations discs



Boehler traction bows, to
be used with STEINMANN
pins, 4mm diam

Boehler traction
bows, to be used
with STEINMANN
pins, 5mm diam

Inside
width

For

Boehler traction bows, to be used with STEINMANN pins, 4mm diam	Boehler traction bows, to be used with STEINMANN pins, 5mm diam	Inside width	For
27-123 with pin (27-135)	27-129 without pin	16x9cm	Arm
27-124 without pin	27-129 without pin		
27-125 with pin (27-137)	27-130 without pin	21x11cm	Femur and knee
27-126 without pin	27-130 without pin		
27-127 with pin (27-139)	27-128 without pin	21x15cm	Femur and knee
27-128 without pin			



27-146
SCHWARZ
Traction bow for fingers,
to be used with KIRSCHNER-
wire up to 1.6 mm diam.



27-147
Small distractor
as a repositioning device for
distraction and compression of small bones,
for wires with 3.0 mm diam.



27-131 - 145
STEINMANN
triangular shaft ended and trocar point

diam	4,0 mm	5,0 mm
Length		
80 mm	27-131	
100 mm	27-132	
120 mm	27-133	
140 mm	27-134	
150 mm	27-135	27-141
160 mm	27-136	
180 mm	27-137	27-142
200 mm	27-138	27-143
210 mm	27-139	
250 mm	27-140	27-144
300 mm		27-145




27-148
Protection cap
for 4,0 mm STEINMANN Pins



27-149
Protection cap
for 5,0 mm STEINMANN Pins







for soft wire up to diam.	0.5 mm	0.6 mm
for hard wire up to diam.	0.3 mm	0.4 mm



for soft wire up to diam.	1.2 mm	1.2 mm
for hard wire up to diam.	0.8 mm	0.8 mm



0.4 mm
0.3 mm



for soft wire up to diam.	1.5 mm
for hard wire up to diam.	1.0 mm



1.5 mm
1.2 mm



1.5 mm
1.0 mm



for soft wire
up to diam. 2,0 mm
For hard wire
up to diam.
specially hardened
cutting edges, front
and lateral cutting
action



for soft wire
up to diam. 2,0 mm
For hard wire
up to diam.
specially hardened
cutting edges, for
cutting off wires close
to the bone



for soft wire
up to diam. 3,0 mm
For hard wire
up to diam. 2,2 mm



for soft wire
up to diam. 2,0 mm
For hard wire
up to diam. 1,6 mm



for soft wire
up to diam. 1,5 mm
For hard wire
up to diam. 1,0 mm



for soft wire
up to diam. 2,0 mm
For hard wire
up to diam. 1,5 mm



for soft wire
up to diam. 3,5 mm
For hard wire
up to diam. 2,5 mm



for soft wires up to 2,5 mm
for hard wires up to 2,0 mm



for soft wires up to 2,5 mm
for hard wires up to 2,0 mm



27-201
320 mm



for soft wire
up to diam.
3.0 mm
For hard wire
up to diam.
3.0 mm



27-202A
25 cm
For wires up to 3.0 mm
Head cutter



27-202B
38 cm
For wires up to 4.0 mm
Head cutter



DUROTIP

27-202
Cutter for cutting off
KIRSCHNER wires,
Tibia head screws,
bolts, pins, up to 6 mm diam.
470 mm



PENNER
27-203A
15.5 cm



27-203
RADOLF
135 mm



27-204
Nail extracting
forceps
145 mm



Drill points only

27-205
Finger nail drill
with 3 drill points



27-209
Nail splitting
forceps
135 mm



27-210
17 cm
Ring cutter



WALTON
27-210A
16 cm
27-210B
Spare blade for ring cutters



27-211
Finger ring saw
150 mm