

SURGICAL

CATALOGUE



PRODUCING QUALITY GOODS FOR 15 YEARS

We are Producers of quality surgical, dental, beauty and general stainless steel products in singleuse and multiuse type. We guarantee on our products and facilitate our customers with easy return policies.



For more information visit www.mjwintl.net



Manufacturers and Exporters of Surgical, Dental & Beauty Instruments

Head Office

BI24S 88/8 Bara Pathar
Christian Town, Sialkot
Pakistan (51310)
0092 524 264829

Branch Office

Flat 1, Kinsale Grange
57 Langley Park Road,
Sutton-Surrey SM2-5GZ UK.
Mob: 00447701084895

Contact Details

Email: info@mjwintl.net
mjwinternational2005@gmail.com
Mobile: 0092 321 712 1510
Mobile: 0092 333 371 3443



FACTORY OVERVIEW



INTRODUCTION

for fifteen years MJW Surgical has been supplying the Healthcare Industry with finest quality surgical instruments & innovative products for surgery. With our extensive knowledge and experience in the medical field, our products are specifically designed to offer both exceptional quality and value for money to Healthcare Professionals and Organisations. During the last few years there have been numerous developments within our company, with the introduction of many new innovative products and services to our portfolio, which has allowed MJW Surgical to make a significant impact on the healthcare market. Our head office is based in Sialkot, Pakistan housing all Customer Service and Administrative functions. Additionally, we have an office in the UK covering Europe and the USA. Our range of specialised surgical instrumentation is uniquely manufactured in Sialkot to the highest standards of craftsmanship. Our goal is to continue to develop new, innovative medical products that will enhance the quality of surgical procedures.

OUR SERVICE COMMITMENT

QUALITY of All surgical instruments and products manufactured and distributed by MJW Surgical conform to the highest European and International Standards with all products bearing the relevant CE mark for their classification.

GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any product proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts.

QUOTATION SERVICE

We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.

PRODUCT INFORMATION

Sectional catalogues are available for most surgical specialities.

INNOVATION

Our objective is to work in close conjunction with all medical and nursing personnel to develop new and innovative products and services.

29-112
For aspirating
fluid from the excretory
ducts of large glands
and for keeping the
operating area clear
170 mm



	A	B	C
29-102	5 mm	1 mm	4.0 mm
29-103	6 mm	1 mm	5.0 mm
29-103A	7 mm	1 mm	5.5 mm
29-103B	8 mm	1 mm	6.0 mm
29-103C	9 mm	1 mm	6.5 mm
29-103D	10 mm	2 mm	7.5 mm
29-103E	12 mm	2 mm	10 mm
29-103F	14 mm	2 mm	11 mm

Uterine suction tube



29-101
13 cm
Uterine suction connector



Uterine suction tubes



29-104
29 cm
Specimen &
aspirating
cannula



29-105
29 cm
Specimen &
aspirating
cannula



Aspirating cannulae
29 cm



29-113
Trocar
For puncturing
150 mm



29-106
8 mm



29-107
8 mm



29-108
10 mm



29-109
12 mm



29-110
14 mm



29-114
For clearing
large quantities
of fluid
190 mm



29-111
12 cm
Coupling piece

29-115 150mm
Suction handle, with vacuum regulator and finger
cut-off, brass, mirror-finish chromium-plated



olive
Ø 9 mm

29-116 170mm
Suction handle, with hand regulation valve and
quick-shutter, brass, mirror-finish chromium-plated



olive
Ø 9 mm

Both suction handles have bayonet fitting
and can be used for the suction cannulae.



29-117
Trocar
For puncturing
170 mm



29-118
Hollow needle,
sharp
150 mm



29-119
double curved
For general
surgery
230 mm




29-120
For tonsillectomy,
oral and laryngeal
operations
190 mm



29-121
For tonsillectomy,
oral and laryngeal
operations
220 mm

With bayonet fitting to be used with
suction handles (29-115 - 29-116)







TUNGESTEN CARBIDE INSTRUMENTS
FOR THOSE.....
WHO DEMAND SECURITY & ACCURACY.



29-210
GENERAL BIOLOGY DISSECTING SET
SINGLE LEATHER CASE.



29-211
BASIC DISSECTING SETS
VINYLE POUCH.



29-212
ANATOMY DISSECTING SET
DOUBLE FOLDING LEATHER CASE.



29-213
LARGE ANIMAL SET
DOUBLE FOLDING
LEATHER CASE.



29-214
MICRO DISSECTING SET
DOUBLE FOLDING
LEATHER CASE.



29-215
NEW INSTRUCTORS DISSECTING SET
(Set Of 12)



29-216
BASIC DISSECTING SETS
In Attractive Zipper Case
(set of 9)



29-217
BASIC DISSECTING SETS (Set of 13)
(In Attractive Zipper Case)



Tungsten Carbide is One Of The Hardest Alloys. We Bond These Inserts Into High Quality German Stainless Steel. Very Special Reasons Of The End Users. We Offer A Wide Range Of Instruments Which Includes Scissors, Forceps And Needle Holders Of All kinds Made Of German Stainless Steel.

Briefly, We Strive To End Up With Quality, Precision And Conformity To Our Buyers/specification And Standard By All Means. The Process Indeed Involves Complex Machining Operations, Skill And Experience.

Let Us Please Secure Your Attention To Serve You Better.

MANAGEMENT