

SURGICAL

CATALOGUE



PRODUCING QUALITY GOODS FOR 15 YEARS

We are Producers of quality surgical, dental, beauty and general stainless steel products in singleuse and multiuse type. We guarantee on our products and facilitate our customers with easy return policies.



For more information visit www.mjwintl.net



Manufacturers and Exporters of Surgical, Dental & Beauty Instruments

Head Office

BI24S 88/8 Bara Pathar
Christian Town, Sialkot
Pakistan (51310)
0092 524 264829

Branch Office

Flat 1, Kinsale Grange
57 Langley Park Road,
Sutton-Surrey SM2-5GZ UK.
Mob: 00447701084895

Contact Details

Email: info@mjwintl.net
mjwinternational2005@gmail.com
Mobile: 0092 321 712 1510
Mobile: 0092 333 371 3443



FACTORY OVERVIEW



INTRODUCTION

for fifteen years MJW Surgical has been supplying the Healthcare Industry with finest quality surgical instruments & innovative products for surgery. With our extensive knowledge and experience in the medical field, our products are specifically designed to offer both exceptional quality and value for money to Healthcare Professionals and Organisations. During the last few years there have been numerous developments within our company, with the introduction of many new innovative products and services to our portfolio, which has allowed MJW Surgical to make a significant impact on the healthcare market. Our head office is based in Sialkot, Pakistan housing all Customer Service and Administrative functions. Additionally, we have an office in the UK covering Europe and the USA. Our range of specialised surgical instrumentation is uniquely manufactured in Sialkot to the highest standards of craftsmanship. Our goal is to continue to develop new, innovative medical products that will enhance the quality of surgical procedures.

OUR SERVICE COMMITMENT

QUALITY of All surgical instruments and products manufactured and distributed by MJW Surgical conform to the highest European and International Standards with all products bearing the relevant CE mark for their classification.

GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any product proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts.

QUOTATION SERVICE

We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.

PRODUCT INFORMATION

Sectional catalogues are available for most surgical specialities.

INNOVATION

Our objective is to work in close conjunction with all medical and nursing personnel to develop new and innovative products and services.

TUNGESTEN CARBIDE INSTRUMENTS
FOR THOSE.....
WHO DEMAND SECURITY & ACCURACY.



29-210
GENERAL BIOLOGY DISSECTING SET
SINGLE LEATHER CASE.



29-211
BASIC DISSECTING SETS
VINYLE POUCH.



29-212
ANATOMY DISSECTING SET
DOUBLE FOLDING LEATHER CASE.



29-213
LARGE ANIMAL SET
DOUBLE FOLDING
LEATHER CASE.



29-214
MICRO DISSECTING SET
DOUBLE FOLDING
LEATHER CASE.



29-215
NEW INSTRUCTORS DISSECTING SET
(Set Of 12)



29-216
BASIC DISSECTING SETS
In Attractive Zipper Case
(set of 9)



29-217
BASIC DISSECTING SETS (Set of 13)
(In Attractive Zipper Case)







Tungsten Carbide is One Of The Hardest Alloys. We Bond These Inserts Into Hight Quality German Stianless SteelVery Special Reasonns Of The End Users. We Offer A Wide Range Of Instruments Which Includes Scissors, Forceps And Needle Holders Of All kinds Made Of German Stainless Steel.
Briefly, We Strive To End Up With Quality, Precision And Conformity To Our Buyers/specification And Standard By All Means. The Process Indeed Involves Complex Machining Operations, Skill And Experience.
Let Us Please Secure Your Atteion To Serve You Better.

MANAGEMENT



Instruments with Tungsten Carbide Inserts

TC GOLD

STANDARD
30-101 to 30-106STANDARD
30-107 to 30-112STANDARD
30-113 to 30-118MAYO
30-119 to 30-124MAYO-LEXER
30-125 to 30-126
19 cmMETZENBAUM
30-127 to 30-133METZENBAUM
30-134 to 30-137IRIS
30-138 to 30-139

Instruments with Tungsten Carbide Inserts

TC GOLD

MAYO-STILLE
30-150 to 30-153METZENBAUM-FINO
30-140 to 30-149LISTER
30-154 to 30-159

30-154 9 cm
30-155 11 cm
30-156 14 cm
30-157 16 cm
30-158 18 cm
30-159 20 cm

BRUNS
30-160
24 cm Serrated

Instruments with Tungsten Carbide Inserts

TC GOLD



Instruments with Tungsten Carbide Inserts

TC GOLD






Instruments with Tungsten Carbide Inserts
TC GOLD



90-456
18 cm / 7"
For soft wires up to 2 mmØ
For hard wires up to 1.5 mmØ



90-457
22 cm / 8 1/2"
For soft wires up to 3.5 mm
For hard wires up to 2.5 mm



TC wire cutter
90-463 38.5 cm / 15" 5 mm
90-464 45 cm / 17 1/2" 6 mm
90-465 58 cm / 22" 7.5 mm



90-458
22 cm / 8 1/2"
For soft wires up to 2.5 mmØ
For hard wires up to 2.2 mmØ



90-459 14 cm / 5 1/2"
90-460 16 cm / 6 1/2"
90-461 18 cm / 7"
90-462 22 cm / 8 1/2"
For soft wires up to 3 mmØ
For hard wires up to 2.2 mmØ

SECTION PAGE	
Surgical Holloware	31
	414 to 424

